

Product datasheet for **AR00100PU-N**

Influenza B (Strain B/Florida/04/06) Protein

Product data:

Product Type:	Native Proteins
Description:	Influenza B Strain B/Florida/04/06 protein, 1 mg
Protein Source:	Chicken eggs
Purity:	>90% pure (SDS-PAGE). Ultracentrifugation with 10-40% sucrose gradient
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 0.05M Tris-HCl, pH 8.0 containing 0.1M Sodium chloride, 5 mM EDTA, 0.09% Sodium azide and 0.005% Thimerosal
Preparation:	Liquid purified protein
Applications:	Serological studies of Influenza B virus, immunogen for antibody production.
Protein Description:	Influenza B Virus, Strain B/Florida/04/06
Note:	Caution: No test guarantees a product to be non-infectious. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product.
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)