

Product datasheet for AR00101PU-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Influenza B Protein

Product data:

Product Type: Native Proteins

Description: Influenza B protein, 1 mg

Protein Source: Chicken eggs

Purity: >90% pure (SDS-PAGE) by Ultracentrifugation with 10-40% sucrose gradient

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 0.05M Tris-HCl, pH 8.0 containing 0.1M Sodium Chloride, 5 mM EDTA, 0.09%

Sodium Azide and 0.005% Thimerosal

Preparation: Liquid purified protein

Applications: Serological studies of Influenza B virus, immunogen for antibody production.

Protein Description: Influenza B Virus, Strain B/Florida/07/04

Note: Caution: No test guarantees a product to be non-infectious. All materials should be handled

as if potentially infectious. Generally accepted laboratory practices appropriate for infectious

materials should be employed when handling this product.

Storage: Store the protein at -20°C.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.