

## Product datasheet for **AR00103PU-N**

### Influenza B (Strain B/Malaysia/2506/04) Protein

#### Product data:

<b>Product Type:</b>	Native Proteins
<b>Description:</b>	Influenza B Strain B/Malaysia/2506/04 protein, 1 mg
<b>Protein Source:</b>	Chicken eggs
<b>Purity:</b>	>90% pure (SDS-PAGE). Ultracentrifugation with 10-40% sucrose gradient
<b>Buffer:</b>	Presentation State: Purified State: Liquid purified protein Buffer System: 0.05M Tris-HCl, pH 8.0 containing 0.1M Sodium chloride, 5 mM EDTA, 0.09% Sodium azide and 0.005% Thimerosal
<b>Preparation:</b>	Liquid purified protein
<b>Applications:</b>	Serological studies of Influenza B virus, immunogen for antibody production.
<b>Protein Description:</b>	Influenza B Virus, Strain B/Malaysia/2506/04
<b>Storage:</b>	Store the protein at -20°C. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.



[View online »](#)