

Product datasheet for **AR00106PU-N**

Influenza B (Strain B/Victoria/504/00) Protein

Product data:

Product Type:	Native Proteins
Description:	Influenza B Strain B/Victoria/504/00 protein, 1 mg
Protein Source:	Chicken eggs
Concentration:	lot specific
Purity:	>90% pure (SDS-PAGE). Ultracentrifugation using 10-40% sucrose gradient
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 0.05M Tris-HCl, pH 8.0 containing 0.1M Sodium chloride, 5 mM EDTA, 0.09% Sodium azide and 0.005% Thimerosal
Preparation:	Liquid purified protein
Applications:	Serological studies of influenza B virus, immunogen for antibody production. Has been tested with anti-influenza B monoclonal antibodies in ELISA.
Protein Description:	Influenza B (B/Victoria/504/00)
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)