

Product datasheet for **AR00134PU-N**

Sm/RNP Antigen Bovine Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Sm/RNP Antigen bovine protein, 0.2 mg
Species:	Bovine
Concentration:	lot specific
Purity:	>90% pure by SDS-PAGE. Purification: Immunoaffinity Chromatography
Buffer:	Presentation State: Aff - Purified State: Liquid purified protein Buffer System: 20 mM HEPES-NaOH, pH 7.5, 0.32 M Potassium Chloride, 10 mM MgCl ₂ , 0.2 mM DTT, 0.1 mM PMSF, 20% Glycerol containing 0.1% preservative, Bromo-Nitro-Dioxane/Methylisothiazoline
Preparation:	Liquid purified protein
Applications:	ELISA: 0.5-1.0 µg/ml coating concentration. After thawing, mix thoroughly by vortexing before use.
Protein Description:	Sm/RNP Antigen
Storage:	Store the protein at -20 or -70°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)