

Product datasheet for **AR00135PU-N**

Sm/RNP Antigen Bovine Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Sm/RNP Antigen bovine protein, 1 kU
Species:	Bovine
Concentration:	lot specific
Purity:	Immobilized human anti-RNP immune globulins. Purification is verified by SDS PAGE. Other antigens (SS-B, SS-A, Scl-70, Jo-1) are not detected with ELISA assays.
Buffer:	Presentation State: Aff - Purified State: Liquid purified Buffer System: 0.01M Tris, 0.53M Sodium chloride containing 0.02% Sodium azide
Bioactivity:	Biological: One unit of activity is the amount of antigen that will generate an absorbance of 1.2 when assayed using standard procedure for enzyme immunoassay. Activity will vary with assay conditions/protocols.
Preparation:	Liquid purified
Applications:	Recommended for use in standard solid phase immunoassays.
Protein Description:	Sm/RNP Complex Antigen
Storage:	Store the protein at -70°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)