

## Product datasheet for **AR00136PU-N**

### SS-A / Ro Bovine Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	SS-A / Ro bovine protein, 0.2 mg
Species:	Bovine
Concentration:	lot specific
Purity:	>90% by SDS-PAGE
Buffer:	Presentation State: Aff - Purified State: Liquid purified Buffer System: 20 mM HEPES-NaOH, pH 7.5, 0.15M Sodium Chloride, 2 mM MgCl <sub>2</sub> , 0.2 mM DTT, 0.1 mM PMSF, 20% Glycerol containing 0.1% preservative and Bromo-Nitro-Dioxane/Methylisothiazoline
Preparation:	Liquid purified
Applications:	ELISA: 0.5-1.0 µg/ml coating concentration. After thawing, mix thoroughly by vortexing before use.
Protein Description:	SS-A (Ro) Antigen
Storage:	Store the protein at -20°C or -70°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_872596</a>
Locus ID:	359715
Summary:	tripartite motif containing 21



[View online »](#)