

## Product datasheet for **AR00137PU-N**

### SS-A / Ro Bovine Protein

#### Product data:

Product Type:	Native Proteins
Description:	SS-A / Ro bovine protein, 1 kU
Species:	Bovine
Protein Source:	Spleen
Concentration:	lot specific
Purity:	Immuno-affinity chromatography. Purification is verified by SDS-PAGE. Other antigens are not detected with ELISA assays.
Buffer:	Presentation State: Aff - Purified State: Liquid purified Buffer System: 0.01M Tris, 0.53M Sodium chloride containing 0.02% Sodium azide
Bioactivity:	Biological: One unit of activity is the amount of antigen that will generate an absorbance of 1.2 when assayed using standard procedure for enzyme immunoassay. Activity will vary with assay conditions/protocols.
Preparation:	Liquid purified
Applications:	Recommended for use in standard solid phase immunoassays.
Protein Description:	SS-A (RO) Antigen
Storage:	Store the protein at -70°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_872596</a>
Locus ID:	359715
Summary:	tripartite motif containing 21



[View online »](#)