

## Product datasheet for **AR00138PU-N**

### SS-B / La Rabbit Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	SS-B / La rabbit protein, 1 kU
Species:	Rabbit
Concentration:	lot specific
Purity:	Purified from rabbit thymus by immuno-affinity chromatography. All rabbit thymus materials are of US origin. Other antigens are not detected with ELISA assays. Purification is verified by SDS PAGE.
Buffer:	Presentation State: Aff - Purified State: Liquid purified Buffer System: 0.01M Tris, 0.53M Sodium chloride containing 0.02% Sodium azide
Bioactivity:	Biological: One unit of activity is the amount of antigen that will generate an absorbance of 1.2 when assayed using standard procedure for enzyme immunoassay. Activity will vary with assay conditions/protocols.
Preparation:	Liquid purified
Applications:	Use in standard solid phase immunoassays.
Protein Description:	Purified SS-B (La) Antigen
Storage:	Store the protein at -70°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	SSB, Sjogren syndrome type B antigen, Lupus La protein, La ribonucleoprotein, La autoantigen, LARP3
Protein Families:	Use in standard solid phase immunoassays.



[View online »](#)