

Product datasheet for **AR01002PU-N**

Mouse IgE Mouse Protein

Product data:

Product Type:	Native Proteins
Description:	Mouse IgE mouse protein, 0.5 mg
Species:	Mouse
Protein Source:	Cell culture
Concentration:	lot specific
Buffer:	Presentation State: Purified State: Liquid purified monoclonal IgE Buffer System: PBS, pH 7.4 containing 0.09% Sodium Azide
Bioactivity:	Specific: Mouse.
Preparation:	Liquid purified monoclonal IgE
Applications:	ELISA: this protein is suitable for use as a Mouse IgE standard in ELISA assays.
Protein Description:	This protein is a purified IgE kappa Mouse monoclonal antibody that recognises Ovalbumin.
Storage:	Store the protein at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)