

## Product datasheet for **AR01003PU-N**

### Mouse IgG1 (kappa chain) Mouse Protein

#### Product data:

|                      |  |
|----------------------|--|
| Product Type:        | Recombinant Proteins   |
| Description:         | Mouse IgG1 (kappa chain) mouse protein, 1 mg   |
| Species:             | Mouse  |
| Concentration:       | lot specific   |
| Purity:              | Single band by agarose electrophoresis and immunoelectrophoresis.  |
| Buffer:              | Presentation State: Purified<br>State: Liquid<br>Buffer System: PBS containing 0.09% Sodium Azide                            |
| Preparation:         | Liquid   |
| Applications:        | <b>ELISA:</b> 0.6-1.25 µg/ml.<br><b>Gel Super Shift Assays.</b>  |
| Protein Description: | Purified Monoclonal Mouse IgG1 kappa (clone MADNP-1) was prepared by affinity chromatography from tissue culture supernatant |
| Storage:             | Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing.         |
| Stability:           | Shelf life: one year from despatch.  |



[View online »](#)