

Product datasheet for **AR01004PU-N**

Mouse IgG2a (kappa chain) Mouse Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Mouse IgG2a (kappa chain) mouse protein, 1 mg
Species:	Mouse
Concentration:	lot specific
Purity:	Affinity chromatography (agarose electrophoresis and immunoelectrophoresis)
Buffer:	Presentation State: Purified State: Liquid Buffer System: PBS containing 0.09% Sodium Azide
Preparation:	Liquid
Applications:	ELISA: 1.25µg/ml.
Protein Description:	Purified mouse IgG2a (Clone: MADNP-2)
Storage:	Store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Mouse Immunoglobulin 2a
Protein Families:	ELISA: 1.25µg/ml.



[View online »](#)