

## Product datasheet for **AR01005PU-N**

### Mouse IgG2b (kappa chain) Mouse Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Mouse IgG2b (kappa chain) mouse protein, 1 mg
Species:	Mouse
Concentration:	lot specific
Purity:	Single band by agarose electrophoresis and immunoelectrophoresis
Buffer:	Presentation State: Purified State: Liquid Buffer System: PBS containing 0.09% Sodium Azide
Preparation:	Liquid
Applications:	ELISA.
Protein Description:	Purified monoclonal IgG2b (Clone: MADNP-3)
Storage:	Store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)