

Product datasheet for **AR03027PU-N**

Hsp90 inhibitor (17-GMB-APA-GA) Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Hsp90 inhibitor (17-GMB-APA-GA) protein, 1 mg
Predicted MW:	767.9
Purity:	>98% 98%
Buffer:	Presentation State: Purified State: Purple solid
Reconstitution Method:	Soluble in DMSO (>25mg/ml) and ethanol (10mg/ml)
Preparation:	Purple solid
Protein Description:	Formula: C ₃₉ H ₅₃ N ₅ O ₁₁
Storage:	Store the antibody (in aliquots) at -20 °C. Can be shipped at 2 - 8 °C. Avoid repeated freezing and thawing. Protect from light!
Stability:	Shelf life: One year from despatch.
Summary:	An analog of geldanamycin containing a maleimido moiety suitable for preparation of geldanamycin immunoconjugates(1-3).



[View online »](#)