

Product datasheet for **AR05013PU-N**

Hepatitis B Core Antigen / HBcAg viral Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Hepatitis B Core Antigen / HBcAg viral recombinant protein, 50 µg
Species:	viral
Concentration:	lot specific
Purity:	>95% by SDS PAGE
Buffer:	State: Liquid purified Buffer System: Phosphate buffered saline containing 50% Glycerol
Preparation:	Liquid purified
Applications:	ELISA.
Protein Description:	This is an <i>E. coli</i> -derived recombinant protein corresponding to full length Hepatitis B core antigen. It is immunoreactive with sera from HBV-infected individuals.
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)