

Product datasheet for **AR05038PU-N**

Low Density Lipoprotein / LDL Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Sample lots of this product are individually evaluated for the labelling of human skin fibroblasts or P-388D cells grown in lipoprotein deficient medium for 48 hours.
Species:	Human
Concentration:	lot specific
Purity:	Purified LDL is labelled with the fluorescent probe Dil (1,1'-dioctadecyl-3,3,3',3'-tetramethyl-indocarbocyanine perchlorate). The resultant product is then dialysed against 0.15 M NaCl, 0.05 M TRIS (pH 7.4) and 0.3 mM EDTA, sterilised by filtration and aseptically packaged.
Preparation:	Liquid purified protein
Storage:	Store the protein undiluted at 2-8°C. DO NOT FREEZE! Shelf life: 3 month from despatch.



[View online »](#)