

Product datasheet for AR05070PU-N

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

Neuropeptide Y / NPY Human Protein

Product data:

Product Type: Recombinant Proteins

Description: Neuropeptide Y / NPY human protein, 0.1 mg

Species: Human

Expression cDNA Clone

or AA Sequence:

Y-P-S-K-P-D-N-P-G-E-D-A-P-A-E-D-M-A-R-Y-Y-S-A-L-R-H-Y-I-N-L-I-T-R-Q-R-Y-NH2

Concentration: lot specific

Purity: >97% pure by HPLC

Buffer: Presentation State: Purified

State: Liquid purified synthetic peptide

Preservative: None Stabilizer: None

Preparation: Liquid purified synthetic peptide

Applications: ELISA.

Protein Description: Synthetic Human Neuropeptide Tyrosine prepared by solid phase synthesis.

Storage: Upon receipt, store undiluted (in aliquots) at -20°C.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: NP 000896

Locus ID: 4852

UniProt ID: <u>P01303</u>, <u>A4D158</u>

Cytogenetics: 7p15.3 Synonyms: PYY4





Neuropeptide Y / NPY Human Protein - AR05070PU-N

Summary:

This gene encodes a neuropeptide that is widely expressed in the central nervous system and influences many physiological processes, including cortical excitability, stress response, food intake, circadian rhythms, and cardiovascular function. The neuropeptide functions through G protein-coupled receptors to inhibit adenylyl cyclase, activate mitogen-activated protein kinase (MAPK), regulate intracellular calcium levels, and activate potassium channels. A polymorphism in this gene resulting in a change of leucine 7 to proline in the signal peptide is associated with elevated cholesterol levels, higher alcohol consumption, and may be a risk factor for various metabolic and cardiovascular diseases. The protein also exhibits antimicrobial activity against bacteria and fungi. [provided by RefSeq, Oct 2014]

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Protein Pathways: Adipocytokine signaling pathway