

Product datasheet for **AR05080PU-N**

Glycogen phosphorylase muscle form Rabbit Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Glycogen phosphorylase muscle form rabbit protein, 25 mg
Species:	Rabbit
Concentration:	lot specific
Purity:	Protein: 90%
Buffer:	Presentation State: Purified State: Liquid purified protein
Bioactivity:	Specific: 60U/mg protein
Preparation:	Liquid purified protein
Applications:	ELISA.
Protein Description:	Native Phosphorylase B from rabbit skeletal muscle. The enzyme preparation contains 5'-AMP which is needed for its catalytic activity.
Note:	1 unit will liberate one μM of glucose-1-phosphate from glycogen and orthophosphate per minute at 37°C and pH 6.8.
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001075653
Locus ID:	100008972
UniProt ID:	P00489
Summary:	phosphorylase, glycogen, muscle



[View online »](#)