

Product datasheet for **AR05089PU-N**

Thyroglobulin Bovine Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Thyroglobulin bovine protein, 0.1 g
Species:	Bovine
Concentration:	lot specific
Buffer:	Presentation State: Purified State: Liquid purified protein
Preparation:	Liquid purified protein
Applications:	ELISA.
Protein Description:	Native Thyroglobulin from bovine thyroid glands is a precursor of the iodated thyroid hormones thyroxine (T4) and triiodothyronine (T3) that affect protein, fat and carbohydrate metabolism. Thyroglobulin is produced by thyrocytes within the thyroid gland and is stored within the follicular lumen.
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_776308
Locus ID:	280706
UniProt ID:	P01266
Cytogenetics:	14q12-q15
Synonyms:	Thyreoglobulin, TGN
Summary:	thyroglobulin
Protein Families:	ELISA.



[View online »](#)