

Product datasheet for **AR05091PU-N**

Trypsin Porcine Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Trypsin porcine protein, 0.5 g
Species:	Porcine
Concentration:	lot specific
Buffer:	State: Partially purified enzyme from porcine pancreas - liquid
Bioactivity:	Specific: 1500 U/mg protein
Reconstitution Method:	For long term storage, addition of a preservative is recommended.
Preparation:	Partially purified enzyme from porcine pancreas - liquid
Applications:	ELISA.
Protein Description:	Native Trypsin from porcine pancreas.
Note:	1 Unit of enzyme causes an increase in absorbance of 0.0003 per minute at 25°C and 253 nm, pH 7.6 on the substrate benzoyl-L arginine ester (BAEE).
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)