

Product datasheet for **AR05098PU-N**

Hepatitis B Core Antigen / HBcAg viral Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Hepatitis B Core Antigen / HBcAg viral recombinant protein, 0.1 mg
Species:	viral
Purity:	~95% by SDS PAGE
Buffer:	State: Liquid purified protein Buffer System: Phosphate buffered saline containing 50% Glycerol
Preparation:	Liquid purified protein
Applications:	ELISA.
Protein Description:	Recombinant Hepatitis B core particles consisting of amino acids 1-183. An <i>Escherichia coli</i> expression system was used.
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)