

Product datasheet for **AR05100PU-N**

Hepatitis B Surface Antigen / HBsAg (ay) viral Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Hepatitis B Surface Antigen / HBsAg (ay) viral recombinant protein, 0.1 mg
Species:	viral
Predicted MW:	24 kDa
Purity:	~95% pure by SDS-PAGE
Buffer:	State: Liquid purified protein. Buffer System: Phosphate buffered saline, pH 7.2 without preservatives
Preparation:	Liquid purified protein.
Applications:	Western Blot ELISA
Protein Description:	Recombinant surface antigen is a AY subtype of the Hepatitis B virus. Genotype A. It is produced using a <i>Saccharomyces cerevisiae</i> . expression system.
Storage:	Store the protein undiluted at 2-8°C. DO NOT FREEZE.
Stability:	Shelf life: one year from despatch.



[View online »](#)