

Product datasheet for **AR05101PU-N**

Hepatitis B Core Antigen / HBcAg viral Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Hepatitis B Core Antigen / HBcAg viral recombinant protein, 0.1 mg
Species:	viral
Concentration:	lot specific
Purity:	95 % by SDS PAGE
Buffer:	State: Liquid purified protein Buffer System: Acetate buffer, pH 4.5
Preparation:	Liquid purified protein
Applications:	Western Blot. ELISA.
Protein Description:	Recombinant Hepatitis B Core produced in <i>Pichia pastoris</i> , amino acids 1-183.
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	HBV Capsid protein, HBV Core protein, p21.5
Protein Families:	Western Blot. ELISA.



[View online »](#)