

## Product datasheet for **AR05103PU-N**

### HCV Core viral Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	HCV Core protein viral recombinant protein, 0.1 mg
Species:	viral
Concentration:	lot specific
Purity:	~95% by SDS PAGE
Buffer:	State: Liquid purified recombinant protein Buffer System: 20 mM Tris 8M Urea; 10 mM betamercaptoethanol
Preparation:	Liquid purified recombinant protein
Applications:	ELISA.
Protein Description:	Recombinant HCV core amino acids 2 to 192 (22 kDa), fused to beta-galactosidase (114 kDa)
Storage:	Store undiluted at -20°C only. Should this product contain a precipitate we recommend microcentrifugation before use.
Stability:	Shelf life: One year from despatch.



[View online »](#)