

Product datasheet for **AR05110PU-N**

Rabbit IgM Rabbit Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Rabbit IgM rabbit protein, 1 mg
Species:	Rabbit
Concentration:	lot specific
Purity:	Purified IgM prepared by a multi-step process including delipidation, selective precipitation and tandem sieve chromatography
Buffer:	Presentation State: Purified State: Liquid purified IgM Buffer System: TRIS buffered saline, pH 8.0 containing 0.09% Sodium Azide (NaN ₃)
Preparation:	Liquid purified IgM
Applications:	ELISA.
Protein Description:	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-rabbit serum and anti-rabbit IgM (μ chain specific). No reaction was observed against anti rabbit IgG F(c). Some light chain cross-reactivity will occur with anti-rabbit IgG.
Note:	Caution: (A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Storage:	Store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)