

## Product datasheet for **AR05132PU-S**

### Plasmodium vivax (circumsporozoite (CSP)) Protozoan Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Plasmodium vivax circumsporozoite (CSP) protozoan recombinant protein, 0.1 mg
Species:	Protozoan
Concentration:	lot specific
Buffer:	State: Liquid purified protein Buffer System: Phosphate buffered saline containing 0.09% Sodium Azide (NaN <sub>3</sub> )
Preparation:	Liquid purified protein
Applications:	ELISA.
Protein Description:	Recombinant <i>Plasmodium vivax</i> CSP protein fused to GST, expressed in <i>E. coli</i> .
Note:	Caution: (A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)