

Product datasheet for **AR05235PU-L**

Bovine IgG1 Bovine Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Bovine IgG1 bovine protein, 1 mg
Species:	Bovine
Concentration:	lot specific
Purity:	Ion Exchange Chromatography
Buffer:	Presentation State: Purified State: Liquid purified IgG1 Buffer System: PBS Preservative: 0.09% Sodium Azide
Preparation:	Liquid purified IgG1
Applications:	ELISA.
Protein Description:	This protein is a preparation of purified Polyclonal Bovine IgG1 prepared from normal Bovine Serum. By Immunoelectrophoresis, it was shown to react with bovine IgG1. No reactivity was observed with antiserum specific for Bovine IgG2, IgA or IgM. It contained < 1% IgG2, IgA and IgM as determined by ELISA.
Note:	Caution: (A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.



[View online »](#)