

Product datasheet for **AR05252PU-N**

Rabbit IgG Rabbit Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Rabbit IgG rabbit protein, 10 mg
Species:	Rabbit
Concentration:	lot specific
Purity:	Ion exchange chromatography from normal rabbit serum
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: Phosphate buffered saline pH 7.2 containing 0.09% Sodium Azide
Preparation:	Liquid purified protein
Applications:	ELISA.
Note:	Caution: (A full Health and Safety assessment is available upon request) This product contains Sodium Azide: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.
Storage:	Store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)