

Product datasheet for **AR08014AG-N**

Bovine IgG (F(ab')₂ Fragment) Bovine Protein

Product data:

Product Type:	Native Proteins
Description:	Bovine IgG F(ab') ₂ Fragment bovine protein, 5 mg
Species:	Bovine
Protein Source:	Serum
Conjugation:	Agarose
Buffer:	State: Liquid serum Buffer System: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 0.01% (w/v) Sodium Azide Presentation Label: Agarose
Preparation:	Liquid serum
Applications:	Suitable for antibody purification or to remove unwanted cross-reactivity.
Protein Description:	Agarose Immobilized Bovine IgG F(ab) ₂ fragment. Sufficient antibody capacity is provided to bind a minimum of 5.0 mg of anti-Bovine IgG.
Storage:	Store the protein undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.



[View online »](#)