

Product datasheet for **AR08015HR-N**

Bovine IgG (Fab Fragment) Bovine Protein

Product data:

Product Type:	Native Proteins
Description:	Bovine IgG Fab Fragment bovine protein, 1 mg
Species:	Bovine
Protein Source:	Serum
Concentration:	lot specific
Purity:	>95% normal serum delipidation, salt fractionation, ion exchange chromatography followed by papain digestion and extensive dialysis against the buffer stated above
Conjugation:	HRP
Buffer:	State: Lyophilized IgG Fab fragment Buffer System: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 10 mg/ml Bovine Serum Albumin (BSA) and 0.01% (w/v) Sodium Azide <u>Label:</u> Horseradish Peroxidase (HRP) Presentation Label: HRP
Reconstitution Method:	Restore with 1.0 ml of deionized water (or equivalent). For extended storage add glycerol to 50%
Preparation:	Lyophilized IgG Fab fragment
Applications:	Western Blot. Dot Blot. ELISA. Immunohistochemistry. Immunoperoxidase electron microscopy.
Protein Description:	Peroxidase Conjugated Bovine IgG Fab fragment.
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)