

Product datasheet for **AR08016FC-N**

Feline IgG (Fc Fragment) Feline Protein

Product data:

Product Type:	Native Proteins
Description:	Feline IgG Fc Fragment feline protein, 1 mg
Species:	Feline
Protein Source:	Serum
Concentration:	lot specific
Conjugation:	FITC
Buffer:	State: Lyophilized IgG F(c) fragment Buffer System: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 10 mg/ml Bovine Serum Albumin (BSA) and 0.01% (w/v) Sodium Azide Label: Fluorescein isothiocyanate (FITC) (Molecular Weight 390 daltons) Protein Ratio: 5.0 moles FITC per mole of Cat IgG F(c) Absorption/Emission: 495 nm / 528 nm Presentation Label: FITC
Reconstitution Method:	Restore with 1.0 ml of deionized water (or equivalent). For extended storage reconstitute product with 50% glycerol instead of water.
Preparation:	Lyophilized IgG F(c) fragment
Applications:	Flow Cytometry. Immunofluorescence.
Protein Description:	Fluorescein Conjugated Cat IgG F(c) fragment.
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the antigen undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)