

Product datasheet for **AR08016HR-N**

Feline IgG (Fc Fragment) Feline Protein

Product data:

Product Type:	Native Proteins
Description:	Feline IgG Fc Fragment feline protein, 1 mg
Species:	Feline
Protein Source:	Serum
Concentration:	lot specific
Conjugation:	HRP
Buffer:	State: Lyophilized IgG F(c) fragment Buffer System: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 10 mg/ml Bovine Serum Albumin (BSA) and 0.01% (w/v) Gentamicin Sulfate Label: Horseradish Peroxidase (HRP) Presentation Label: HRP
Reconstitution Method:	Restore with 1.0 ml of deionized water (or equivalent). For extended storage, mix with glycerol to 50%.
Preparation:	Lyophilized IgG F(c) fragment
Applications:	Western Blot. ELISA. Immunohistochemistry. Immunoperoxidase electron microscopy.
Protein Description:	Peroxidase Conjugated Cat IgG F(c) fragment.
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the antigen undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)