

Product datasheet for **AR08018HR-N**

Feline IgG (Fab fragment) Feline Protein

Product data:

| | |
|------------------------|---|
| Product Type: | Native Proteins |
| Description: | Feline IgG Fab fragment feline protein, 1 mg |
| Species: | Feline |
| Protein Source: | Serum |
| Concentration: | lot specific |
| Conjugation: | HRP |
| Buffer: | State: Lyophilized IgG Fab fragment Buffer System: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 10 mg/ml Bovine Serum Albumin (BSA) and 0.01% (w/v) Gentamicin Sulfate Label: Horseradish Peroxidase (HRP) Presentation Label: HRP |
| Reconstitution Method: | Restore with 1.0 ml of deionized water (or equivalent). For extended storage mix with glycerol to 50%. |
| Preparation: | Lyophilized IgG Fab fragment |
| Applications: | Western Blot. ELISA. Immunohistochemistry. Immunoperoxidase electron microscopy. |
| Protein Description: | Peroxidase Conjugated Cat IgG Fab fragment |
| Storage: | Prior to reconstitution store at 2-8°C. Following reconstitution store the antigen undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |



[View online »](#)