

## Product datasheet for **AR08018TC-N**

### Feline IgG (Fab fragment) Feline Protein

#### Product data:

Product Type:	Native Proteins
Description:	Feline IgG Fab fragment feline protein, 1 mg
Species:	Feline
Protein Source:	Serum
Concentration:	lot specific
Conjugation:	Rhodamine
Buffer:	State: Lyophilized IgG Fab fragment Buffer System: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 10 mg/ml Bovine Serum Albumin (BSA) and 0.01% (w/v) Sodium Azide <b>Label:</b> Tetramethylrhodamine isothiocyanate (TRITC) (Molecular Weight 444 daltons) <b>Protein Ratio:</b> 1.2 moles TRITC per mole of Cat IgG Fab <b>Absorption/Emission</b> 550 nm / 570 nm Presentation Label: Rhodamine
Reconstitution Method:	Restore with 1.0 ml of deionized water (or equivalent). For extended storage reconstitute product with 50% glycerol instead of water.
Preparation:	Lyophilized IgG Fab fragment
Applications:	Flow Cytometry. Immunofluorescence.
Protein Description:	Rhodamine Conjugated Cat IgG Fab fragment
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the antigen undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)