

Product datasheet for **AR08020AG-N**

Chicken IgG / Chicken IgY (F(ab')₂ fragment) Chicken Protein

Product data:

Product Type:	Native Proteins
Description:	Chicken IgG / Chicken IgY F(ab') ₂ fragment chicken protein, 5 mg
Species:	Chicken
Protein Source:	Serum
Purity:	>95% by chromatography coupled to activated agarose
Conjugation:	Agarose
Buffer:	State: Liquid IgG F(ab) ₂ fragment Buffer System: 0.01 M Sodium Phosphate, 0.14 M Sodium Chloride, pH 7.4 containing 0.09% (w/v) Sodium Azide Presentation Label: Agarose
Preparation:	Liquid IgG F(ab) ₂ fragment
Applications:	Suitable for antibody purification or to remove unwanted cross-reactivity.
Protein Description:	Agarose Immobilized Chicken IgG F(ab) ₂ fragment
Note:	<u>Physical State</u> : Suspension of agarose beads <u>Antibody/Agarose</u> : 1 mg antibody per cc agarose
Storage:	Store the antigen undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.



[View online »](#)