

Product datasheet for **AR08022PU-N**

Chicken IgM (Mu chain) Chicken Protein

Product data:

Product Type:	Native Proteins
Description:	Chicken IgM (Mu chain) chicken protein, 1 mg
Species:	Chicken
Protein Source:	Serum
Concentration:	lot specific
Purity:	>95% prepared from normal serum by a multi-step process which includes delipidation, selective precipitation and tandem molecular sieve chromatography followed by extensive dialysis against the buffer stated above
Buffer:	Presentation State: Purified State: Liquid Buffer System: 0.1 M Tris Chloride, 0.5 M Sodium Chloride, pH 8.0 containing 0.1% (w/v) Sodium Azide
Preparation:	Liquid
Applications:	Suitable for use as antigen or ligand in immunochemical reactions, as a control or standard in assays, for conjugation and most other immunological methods requiring highly purified immunoglobulins.
Protein Description:	Chicken IgM
Storage:	Store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)