

Product datasheet for **AR09082PU-N**

RHEB (1-181, T7-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	RHEB (1-181, T7-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MASMTGGQQM GRGSASMPQS KSRKIAILGY RSVGKSSLTI QFVEGQFVDS YDPTIENTFT KLITVNGQEY HLQLVDTAGQ DEYSIFPQTY SIDINGYILV YSVTSIKSFE VIKVIHGKLL DMVGKQVQPI MLVGNKKDLH MERVISYEEG KALAESWNAA FLESAKENQ TAVDVFRRRI LEAEKMDGAA SQGKSSC
Tag:	T7-tag
Predicted MW:	21.7 kDa
Concentration:	lot specific
Purity:	>95% by SDS-PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris pH 8.0, 1 mM DTT, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human RheB protein, fused to T7-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_005605
Locus ID:	6009
UniProt ID:	Q15382 , A0A090N900
Cytogenetics:	7q36.1
Synonyms:	RHEB2



[View online »](#)

Summary:

This gene is a member of the small GTPase superfamily and encodes a lipid-anchored, cell membrane protein with five repeats of the RAS-related GTP-binding region. This protein is vital in regulation of growth and cell cycle progression due to its role in the insulin/TOR/S6K signaling pathway. The protein has GTPase activity and shuttles between a GDP-bound form and a GTP-bound form, and farnesylation of the protein is required for this activity. Three pseudogenes have been mapped, two on chromosome 10 and one on chromosome 22. [provided by RefSeq, Jul 2008]

Protein Pathways:

Insulin signaling pathway, mTOR signaling pathway

Product images: