

Product datasheet for **AR09130PU-L**

BMP7 (293-431, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	BMP7 (293-431, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Tag:	His-tag
Predicted MW:	16.8 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 10 mM Sodium Citrate buffer (pH 3.5) containing 10% glycerol
Endotoxin:	< 1.0 EU per 1 µg of protein (determined by LAL method)
Preparation:	Liquid purified protein
Protein Description:	Recombinant human BMP-7 protein, fused to His-tag, was expressed in E.coli and purified by using conventional chromatography techniques, after refolding of the isolated inclusion body in renaturation buffer.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001710
Locus ID:	655
Cytogenetics:	20q13.31
Synonyms:	OP-1



[View online »](#)

Summary:

This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer, which plays a role in bone, kidney and brown adipose tissue development. Additionally, this protein induces ectopic bone formation and may promote fracture healing in human patients. [provided by RefSeq, Jul 2016]

Protein Families:

Adult stem cells, Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Secreted Protein, Stem cell relevant signaling - TGFb/BMP signaling pathway

Protein Pathways:

Cytokine-cytokine receptor interaction, Hedgehog signaling pathway, TGF-beta signaling pathway

Product images: