

Product datasheet for **AR09167PU-L**

SUMO2 (1-93) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	SUMO2 (1-93) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MADEKPKEGV K TENNDHINL KVAGQDGSVW QFKIKRHTPL SKLMKAYCER QGLSMRQIRF RFDGQPINET DTPAQLEMED EDTIDVFQQQ TGG
Predicted MW:	10.6 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified peptide Buffer System: 20 mM Tris-HCl buffer (pH 8.0)
Preparation:	Liquid purified peptide
Protein Description:	Recombinant Human SUMO2 protein was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001005849
Locus ID:	6613
UniProt ID:	P61956
Cytogenetics:	17q25.1
Synonyms:	HSMT3; Smt3A; SMT3B; SMT3H2; SUMO3



[View online »](#)

Summary:

This gene encodes a protein that is a member of the SUMO (small ubiquitin-like modifier) protein family. It functions in a manner similar to ubiquitin in that it is bound to target proteins as part of a post-translational modification system. However, unlike ubiquitin which targets proteins for degradation, this protein is involved in a variety of cellular processes, such as nuclear transport, transcriptional regulation, apoptosis, and protein stability. It is not active until the last two amino acids of the carboxy-terminus have been cleaved off. Numerous pseudogenes have been reported for this gene. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome

Product images: