

## Product datasheet for AR09182PU-L

## 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

## FABP3 (1-133, His-tag) Human Protein

**Product data:** 

**Product Type: Recombinant Proteins** 

**Description:** FABP3 (1-133, His-tag) human recombinant protein, 0.5 mg

Species: Human E. coli **Expression Host:** 

**Expression cDNA Clone** 

MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSHMVD AFLGTWKLVD SKNFDDYMKS or AA Sequence: LGVGFATRQV ASMTKPTTII EKNGDILTLK THSTFKNTEI SFKLGVEFDE TTADDRKVKS IVTLDGGKLV

HLQKWDGQET TLVRELIDGK LILTLTHGTA VCTRTYEKEA

Tag: His-tag Predicted MW: 19.1 kDa Concentration: lot specific

**Purity:** >95% by SDS - PAGE

**Buffer:** Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

**Endotoxin:** < 1.0 EU per 1 µg of protein (determined by LAL method )

Liquid purified protein Preparation:

**Protein Description:** Recombinant FABP3 protein was expressed in E.coli and purified by using conventional

chromatography techniques.

Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Storage:

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

NP 001307925 RefSeq:

Locus ID: 2170 **Cytogenetics:** 1p35.2

FABP11; H-FABP; M-FABP; MDGI; O-FABP Synonyms:





**Summary:** 

The intracellular fatty acid-binding proteins (FABPs) belongs to a multigene family. FABPs are divided into at least three distinct types, namely the hepatic-, intestinal- and cardiac-type. They form 14-15 kDa proteins and are thought to participate in the uptake, intracellular metabolism and/or transport of long-chain fatty acids. They may also be responsible in the modulation of cell growth and proliferation. Fatty acid-binding protein 3 gene contains four exons and its function is to arrest growth of mammary epithelial cells. This gene is a candidate tumor suppressor gene for human breast cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016]

**Protein Pathways:** 

PPAR signaling pathway

## **Product images:**

