

## Product datasheet for **AR09241PU-L**

### VAMP-4 (1-115, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	VAMP-4 (1-115, His-tag) human recombinant protein, 0.25 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MPPKFKRHLN DDDVTGSVKS ERRNLEDDS DEEEDFFLRG PSGPRFGPRN DKIKHVQNQV DEVIDVMQEN ITKVIERGER LDELQDKSES LSDNATAFSN RSKQLRRQMW WRGCKLEHHH HHH
Tag:	His-tag
Predicted MW:	14.5 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.1 M NaCl, and 20% glycerol
Endotoxin:	< 1.0 EU per 1 µg of protein (determined by LAL method )
Preparation:	Liquid purified protein
Protein Description:	Recombinant human VAMP4 protein, fused to His-tag at C-terminus, was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_001172056</a>
Locus ID:	8674
UniProt ID:	<a href="#">O75379</a>
Cytogenetics:	1q24.3
Synonyms:	VAMP-4; VAMP24



[View online »](#)

**Summary:**

Synaptobrevins/VAMPs, syntaxins, and the 25-kD synaptosomal-associated protein SNAP25 are the main components of a protein complex involved in the docking and/or fusion of synaptic vesicles with the presynaptic membrane. The protein encoded by this gene is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family. This protein may play a role in trans-Golgi network-to-endosome transport. [provided by RefSeq, Jul 2008]

**Protein Families:**

Transmembrane

**Protein Pathways:**

SNARE interactions in vesicular transport

**Product images:**