

Product datasheet for AR09245PU-N

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

SAT1 (1-171, His-tag) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: SAT1 (1-171, His-tag) human recombinant protein, 0.1 mg

Species: Human E. coli **Expression Host:**

Expression cDNA Clone

MGSSHHHHHH SSGLVPRGSH MAKFVIRPAT AADCSDILRL IKELAKYEYM EEQVILTEKD or AA Sequence:

LLEDGFGEHP FYHCLVAEVP KEHWTPEGHS IVGFAMYYFT YDPWIGKLLY LEDFFVMSDY RGFGIGSEIL KNLSQVAMRC RCSSMHFLVA EWNEPSINFY KRRGASDLSS EEGWRLFKID

KEYLLKMATE E

Tag: His-tag Predicted MW: 22.1 kDa **Concentration:** lot specific

>95% by SDS - PAGE **Purity:**

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Preparation: Liquid purified protein

Protein Description: Recombinant human SAT-1, fused to His-tag at N-terminus, was expressed in E.coli and

purified by using conventional chromatography.

Storage: Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer.

Avoid repeated freezing and thawing.

Shelf life: one year from despatch. Stability:

RefSeq: NP 002961

6303 Locus ID:

UniProt ID: P21673, A0A384NQ10

Cytogenetics: Xp22.11

Synonyms: DC21; KFSD; KFSDX; SAT; SSAT; SSAT-1





Summary: The protein encoded by this gene belongs to the acetyltransferase family, and is a rate-

limiting enzyme in the catabolic pathway of polyamine metabolism. It catalyzes the

acetylation of spermidine and spermine, and is involved in the regulation of the intracellular concentration of polyamines and their transport out of cells. Defects in this gene are

associated with keratosis follicularis spinulosa decalvans (KFSD). Alternatively spliced

transcripts have been found for this gene.[provided by RefSeq, Sep 2009]

Protein Families: Druggable Genome

Protein Pathways: Arginine and proline metabolism, Metabolic pathways

Product images:

