

Product datasheet for **AR09262PU-L**

EIF4EBP1 / 4E-BP1 (1-118, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	EIF4EBP1 / 4E-BP1 (1-118, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSGLVPRGSH</u> MSGSSCSQT PSRAIPATRR VVLGDGVQLP PGDYSTTPGG TLFSTTPGGT RIIYDRKFLM ECRNSPVTKT PPRDLPTIPG VTSPSSDEPP MEASQSHLRN SPEDKRAGGE ESQFEMDI
Tag:	His-tag
Predicted MW:	14.7 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: PBS (pH 7.4) containing 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant EIF4EBP1 protein, fused to His-tag, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_004086</u>
Locus ID:	1978
UniProt ID:	<u>Q13541</u>
Cytogenetics:	8p11.23
Synonyms:	4E-BP1; 4EBP1; BP-1; PHAS-I



[View online »](#)

Summary:

This gene encodes one member of a family of translation repressor proteins. The protein directly interacts with eukaryotic translation initiation factor 4E (eIF4E), which is a limiting component of the multisubunit complex that recruits 40S ribosomal subunits to the 5' end of mRNAs. Interaction of this protein with eIF4E inhibits complex assembly and represses translation. This protein is phosphorylated in response to various signals including UV irradiation and insulin signaling, resulting in its dissociation from eIF4E and activation of mRNA translation. [provided by RefSeq, Jul 2008]

Protein Pathways:

Acute myeloid leukemia, ErbB signaling pathway, Insulin signaling pathway, mTOR signaling pathway

Product images: