

Product datasheet for **AR09304PU-N**

SMAD3 (1-425, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	SMAD3 (1-425, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Tag:	His-tag
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl and 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant SMAD3 protein, fused to His-tag, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001138574
Locus ID:	4088
Cytogenetics:	15q22.33
Synonyms:	HSPC193; HsT17436; JV15-2; LDS1C; LDS3; MADH3
Summary:	The SMAD family of proteins are a group of intracellular signal transducer proteins similar to the gene products of the Drosophila gene 'mothers against decapentaplegic' (Mad) and the C. elegans gene Sma. The SMAD3 protein functions in the transforming growth factor-beta signaling pathway, and transmits signals from the cell surface to the nucleus, regulating gene activity and cell proliferation. It also functions as a tumor suppressor. Mutations in this gene are associated with aneurysms-osteoarthritis syndrome and Loeyes-Dietz Syndrome 3. [provided by RefSeq, Nov 2019]



[View online »](#)

Protein Families: Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Stem cell relevant signaling - JAK/STAT signaling pathway, Stem cell relevant signaling - TGFb/BMP signaling pathway, Transcription Factors

Protein Pathways: Adherens junction, Cell cycle, Chronic myeloid leukemia, Colorectal cancer, Pancreatic cancer, Pathways in cancer, TGF-beta signaling pathway, Wnt signaling pathway

Product images:

