

Product datasheet for **AR09308PU-L**

MAT2B (1-323) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	MAT2B (1-323) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MPEMPEDMEQ EEVNIPNRRV LVTGATGLLG RAVHKEFQQN NWHAVGCGFR RARPKFEQVN LLDSNAVHHI IHDFQPHVIV HCAAERRPDV VENQPDAASQ LNVDAAGNLA KEAAAVGAFL IYSSDYVFD GTNPPYREED IPAPLNLYGK TKLDGEKAVL ENNLGAAVLR IPILYGEVEK LEESAIVTMF DKVQFSNKSA NMDHWQQRFP THVKDVATVC RQLAEKRMLD PSIKGTFHWS GNEQMTKYEM ACAIADAFNL PSSHLRPITD SPVLGAQRPR NTQLDCSKLE TLGIGQRTPF RIGIKESLWP FLIDKRWRQT VFH
Predicted MW:	36.4 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 1 mM DTT, 1 mM EDTA, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human MAT II beta was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_037415
Locus ID:	27430
UniProt ID:	Q9NZL9 , A0A140VJP2
Cytogenetics:	5q34
Synonyms:	MAT-II; MATIIbeta; Nbla02999; SDR23E1; TGR



[View online »](#)

Summary:

The protein encoded by this gene belongs to the methionine adenosyltransferase (MAT) family. MAT catalyzes the biosynthesis of S-adenosylmethionine from methionine and ATP. This protein is the regulatory beta subunit of MAT. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Nov 2012]

Protein Pathways:

Cysteine and methionine metabolism, Metabolic pathways, Selenoamino acid metabolism

Product images: