

## Product datasheet for **AR09451PU-N**

### Protein CutA (33-179, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Protein CutA (33-179, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MRLLLLPRVL LTMASGSPPT QPSPASDSGS GYVPGSVSAA FVTCPNKVA KEIARAVVEK RLAACVNLIQ QITSIYEWKG KIEEDSEVLM MIKTQSSLVP ALTDFVRSVH PYEVAEVIAl PVEQGNFPYL QWVRQVTESV SDSITVLP <u>LE HHHHHH</u>
Tag:	His-tag
Predicted MW:	17.1 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 1 mM DTT, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant Human CUTA protein, fused to <i>His-tag</i> at C-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u><a href="#">NP_001014433</a></u>
Locus ID:	51596
UniProt ID:	<u><a href="#">O60888</a></u>
Cytogenetics:	6p21.32
Synonyms:	ACHAP; C6orf82
Summary:	May form part of a complex of membrane proteins attached to acetylcholinesterase (AChE). [UniProtKB/Swiss-Prot Function]



[View online »](#)

Protein Families: Transmembrane

**Product images:**

