

Product datasheet for **AR09473PU-L**

FHIT (1-147, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	FHIT (1-147, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MSFRFGQHLLI KPSVWFLKTE LSFALVNRKP VVPGHVLVCP LRPVERFHDLL RPDEVADLFLQ TTQRVGTVVE KHFGTSLTF SMQDGPEAGQ TVKHVHVHVL PRKAGDFHRN DSIYEELQKH DKEDFPASWR SEEEMAAEAA ALRVYFQLEH <u>HHHHH</u>
Tag:	His-tag
Predicted MW:	17.9 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl Buffer (pH 8.0) containing 10% Glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human FHIT, fused to His-tag at C-terminus, was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_001159715</u>
Locus ID:	2272
UniProt ID:	<u>P49789</u> , <u>A0A024R366</u>
Cytogenetics:	3p14.2
Synonyms:	AP3Aase; FRA3B



[View online »](#)

Summary:

The protein encoded by this gene is a P1-P3-bis(5'-adenosyl) triphosphate hydrolase involved in purine metabolism. This gene encompasses the common fragile site FRA3B on chromosome 3, where carcinogen-induced damage can lead to translocations and aberrant transcripts. In fact, aberrant transcripts from this gene have been found in about half of all esophageal, stomach, and colon carcinomas. The encoded protein is also a tumor suppressor, as loss of its activity results in replication stress and DNA damage. [provided by RefSeq, Aug 2017]

Protein Pathways:

Non-small cell lung cancer, Purine metabolism, Small cell lung cancer

Product images: