

Product datasheet for **AR09526PU-N**

skp (21-161, His-tag) Escherichia coli Protein


Product data:

Product Type:	Recombinant Proteins
Description:	skp (21-161, His-tag) e. coli recombinant protein, 0.1 mg
Species:	Escherichia coli
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SGLVPRGSH</u> MADKIAIVNM GSLFQQVAQK TGVSNLENE FKGRASELQR METDLQAKMK KLQSMKAGSD RTKLEKDVMA QRQTFAQKAQ AFEQDRARRS NEERGLVTR IQTAVKSVAN SQDIDLVDVA NAVAYNSSDV KBITADVLKQ VK
Tag:	His-tag
Predicted MW:	17.9 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 20% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant E.coli Skp, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Summary:	The Seventeen Kilodalton Protein (Skp) is a trimeric periplasmic chaperone that assists outer membrane proteins in their folding and insertion into membranes. This protein required for the efficient release of ompA from the inner membrane, the maintenance of its solubility in the periplasm, and, in association with lipopolysaccharide (LPS), for the efficient folding and insertion of ompA into the outer membrane.



[View online »](#)

Product images:



15% SDS-PAGE (3ug)