

## Product datasheet for **AR09555PU-N**

### L-xylulose reductase (1-244, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	L-xylulose reductase (1-244, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSGLVPRGSH</u> MELFLAGRRV LVTGAGKIG RGTVQALHAT GARVAVSRT QADLDSLVRE CPGIEPVCVD LGDWEATERA LGSVGPVDLL VNNAVALLQ PFLEVTKEAF DRSEFVNLR VQVSQIVAR GLIARGVPGA IVNVSSQCSQ RAVTNHSVYC STKGALDMLT KVMALELGPH KIRVNAVNP VVMTSMGQAT WSDPHKAKTM LNRIPLGKFA EVEHVVNAIL FLLSDRSGMT TGSTLPVEGG FWAC
Tag:	His-tag
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 1 mM DTT, 20% glycerol, 50 mM NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant human DCXR, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_001182147</u>
Locus ID:	51181
UniProt ID:	<u>Q7Z4W1</u>
Cytogenetics:	17q25.3
Synonyms:	DCR; HCR2; HCR2; KIDCR; P34H; PNTSU; SDR20C1; XR



[View online »](#)

**Summary:**

The protein encoded by this gene acts as a homotetramer to catalyze diacetyl reductase and L-xylulose reductase reactions. The encoded protein may play a role in the uronate cycle of glucose metabolism and in the cellular osmoregulation in the proximal renal tubules. Defects in this gene are a cause of pentosuria. Two transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Aug 2010]

**Protein Families:**

Druggable Genome

**Protein Pathways:**

Metabolic pathways, Pentose and glucuronate interconversions

**Product images:**