

## Product datasheet for **AR09660PU-N**

### Secretagogin (SCGN) (1-276, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Secretagogin (SCGN) (1-276, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSGLVPRGSH</u> MDSSREPTLG RLDAAGFWQV WQRFDADEKG YIEEKELDAF FLHMLMKLGT DDTVMKANLH KVKQQFMTTQ DASKDGRIRM KELAGMFLSE DENFLLFRF ENPLDSSVEF MQIWRKYDAD SSGFISAAEL RNFLRDLFLH HKKAISEAKL EEYTGTMKI FDRNKDGRLD LNDLARILAL QENFLQFKM DACSTEERKR DFEKIFAYYD VSKTGALEGP EVDGFVKDMM ELVQPSISGV DLDKFREILL RHCDVKNKDGK IQKSELALCL GLKINP
Tag:	His-tag
Predicted MW:	34.2 kDa
Concentration:	lot specific
Purity:	>90% by SDS-PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1 mM DTT, 0.1M NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant SCGN protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u><a href="#">NP_008929</a></u>
Locus ID:	10590
UniProt ID:	<u><a href="#">O76038</a></u>
Cytogenetics:	6p22.2
Synonyms:	CALBL; DJ501N12.8; SECRET; SEGN; setagin



[View online »](#)

**Summary:**

The encoded protein is a secreted calcium-binding protein which is found in the cytoplasm. It is related to calbindin D-28K and calretinin. This protein is thought to be involved in KCL-stimulated calcium flux and cell proliferation. [provided by RefSeq, Jul 2008]

**Product images:**